




CNTFR Antibody

Product Code	CSB-PA276708
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26992
Immunogen	Peptide sequence around aa.119~123(N-T-Y-P-K) derived from Human CNTFR α
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CNTFR α protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MGC1774
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNTFR

Image



Western blot analysis of extracts from Rat brain tissue and HeLa cells using CNTFR α Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.