




KRT7 Antibody

Product Code	CSB-PA128660
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08729
Immunogen	Peptide sequence around aa.22~26(A-Q-V-R-L) derived from Human Keratin 7
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Keratin 7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	K7; CK7; SCL; K2C7; MGC3625
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KRT7

Image



Western blot analysis of extracts from HeLa cells using Keratin 7 Antibody.