




RPS6KA5 (Ab-376) Antibody

Product Code	CSB-PA144527
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75582
Immunogen	Peptide sequence around aa.374~378 (G-Y-S-F-V) derived from Human MSK1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total MSK1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MSPK1; RPS6KA5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA5

Image



Western blot analysis of extracts from HUVEC cells using MSK1(Ab-376) Antibody and the same antibody preincubated with blocking peptide.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.