





SHC1 (Ab-349) Antibody

Product Code	CSB-PA208610
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29353
Immunogen	Peptide sequence around aa.347~351 (H-Q-Y-Y-N) derived from Human Shc1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total Shc1 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	SH2 domain protein C1; SHC; SHC-transforming protein 1; SHCA; Src homology 2 domain-containing-transforming protein C1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SHC1

Image



Western blot analysis of extracts from HepG2 cells using Shc1 (Ab-349) Antibody.



Immunofluorescence staining of methanol-fixed Hela cells using Shc1 (Ab-349) Antibody.

