




CASP8 Antibody

Product Code	CSB-PA982156
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14790
Immunogen	Peptide sequence around aa.217~221 (S-E-S-Q-T) derived from Caspase 8.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total Caspase 8 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	ICE-like apoptotic protease 5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CASP8

Image



Western blot analysis of extract from 293, HeLa and HT29 cells using Caspase8 Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.