




HNF1B Antibody

Product Code	CSB-PA151885
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35680
Immunogen	Peptide sequence around aa.252~256(R-Q-K-N-P) derived from Human HNF-1b(TCF-2).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total HNF-1b(TCF-2) protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	FJHN; HNF2; LFB3; TCF2; HPC11
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HNF1B

Image



Western blot analysis of extracts from HeLa, K562 and 3T3 cells using HNF-1b(TCF-2) Antibody.