




IGFBP3 Antibody

Product Code	CSB-PA799945
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17936
Immunogen	Peptide sequence around aa.185~189(R-Y-K-V-D) derived from Human IGFBP-3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total IGFBP-3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	IBP3; BP-53;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IGFBP3

Image



Western blot analysis of extracts from various cells using IGFBP-3 Antibody.