




IGFBP3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA799945 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P17936 |
| Immunogen | Peptide sequence around aa.185~189(R-Y-K-V-D) derived from Human IGFBP-3. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous level of total IGFBP-3 protein. |
| Tested Applications | ELISA, WB; WB: 1:500-1:1000 |
| Form | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide. |
| Clonality | Polyclonal |
| Alias | IBP3; BP-53; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IGFBP3 |

Image



Western blot analysis of extracts from various cells using IGFBP-3 Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.