




LGR5 Antibody

Product Code	CSB-PA622827
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75473
Immunogen	Peptide sequence around aa.490~494(K-I-S-N-Q) derived from Human LGR5(GPR49).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total LGR5(GPR49) protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	FEX; GPR49; GPR67; GRP49; MGC117008
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LGR5

Image



Western blot analysis of extracts from HeLa, 3T3 and HepG2 cells using LGR5(GPR49) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.