




SS18 Antibody

Product Code	CSB-PA173100
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15532
Immunogen	Peptide sequence around aa.14~18(G-E-I-T-P) derived from Human SYT.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total SYT protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	SSXT; MGC116875;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SS18

Image



Western blot analysis of extracts from A431 cells using SYT Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.