




LCK (Ab-394) Antibody

Product Code	CSB-PA977883
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06239
Immunogen	Peptide sequence around aa.392-396(N-E-Y-T-A) derived from Human Lck.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Lck protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	LSK; YT16; p56lck; pp58lck;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LCK

Image



Western blot analysis of extracts from JK cells using Lck(Ab-394) Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.