




# ACTB Antibody

<b>Product Code</b>	CSB-PA551635
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P60709
<b>Immunogen</b>	Peptide sequence around aa.6-10(A-A-L-V-V)derived from Human beta-actin.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:1000
<b>Clonality</b>	Polyclonal
<b>Alias</b>	BRWS1, PS1TP5BP1,
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ACTB

## Image



Western blot analysis of extracts from 293, HeLa and 3T3 cells using  $\beta$ -actin Antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.