




CAV2 Antibody

Product Code	CSB-PA901158
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51636
Immunogen	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total Caveolin-2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Clonality	Polyclonal
Alias	CAV,
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CAV2

Image



Western blot analysis of extracts from HeLa and HUVEC cells using Caveolin-2 Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.