




SMAD1 (Ab-206) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA108987 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q15797 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB;WB:1:500-1:1000 |
| Clonality | Polyclonal |
| Alias | BSP1, JV41, MADH1 , MADR1, |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SMAD1 |

Image



Western blot analysis of extracts from HeLa and 3T3 cells using Smad1 (Ab-206) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.