




# Phospho-ACACA (Ser79) Antibody

<b>Product Code</b>	CSB-PA937691
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q13085
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	;
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ACC, ACAC, ACC1, ACCA,
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ACACA

## Image



Western blot analysis of extracts from 293 cells untreated or treated with AICAR using Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody.

<b>Product Modify</b>	Phospho-Ser79
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.