





Phospho-PTEN (Ser380) Antibody

Product Code	CSB-PA585981
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60484
Immunogen	Peptide sequence around phosphorylation site of serine 380(R-Y-S(p)-D-T) derived from Human PTEN.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC, IF; IHC:1:50-1:100, IF:1:100-1:200
Form	Liquid
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTEN

Image



Immunofluorescence staining of methanol-fixed HeLa cells using PTEN (phospho-Ser380) antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using PTEN (phospho-Ser380) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.