





# Phospho-ABL1 (Tyr412) Antibody

<b>Product Code</b>	CSB-PA980598
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P00519
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	ELISA, IHC, IF; IHC:1:50-1:100, IF:1:100-1:200
<b>Form</b>	Liquid
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ABL, JTK7, p150, c-ABL,
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ABL1

## Image



Immunofluorescence staining of methanol-fixed HeLa cells using c-Abl (phospho-Tyr412) antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-Abl (phospho-Tyr412) antibody.

**Product Modify** Phospho-Tyr412