





Phospho-ABL1 (Tyr412) Antibody

Product Code	CSB-PA980598
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00519
Immunogen	Synthesized phosho peptide around human ABL1 412 (L-A-Y(p)-N-K)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC, IF; IHC:1:50-1:100, IF:1:100-1:200
Form	Liquid
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABL1

Image



Immunofluorescence staining of methanol-fixed HeLa cells using c-Abl (phospho-Tyr412) antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-Abl (phospho-Tyr412) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.