




Phospho-MET (Tyr1234/1235) Antibody

Product Code	CSB-PA153714
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08581
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	;
Clonality	Polyclonal
Alias	HGFR, AUTS9, RCCP2, c-Met ,
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MET

Image



Western blot analysis of extracts from HUVEC cells untreated or treated with EGF using Met (Phospho-Tyr1234/1235) Antibody.

Product Modify	Phospho-Tyr1234/1235
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.