




BBC3 Antibody

Product Code	CSB-PA211852
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXH1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; WB:1:500-1:1000
Form	Liquid
Clonality	Polyclonal
Alias	BBC3, JFY1, JFY-1,
Immunogen Species	Homo sapiens (Human)
Target Names	BBC3

Image



Western blot analysis of extracts from 293 and MCF-7 cells using Puma Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.