




# BBC3 Antibody

<b>Product Code</b>	CSB-PA211852
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BXH1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; WB:1:500-1:1000
<b>Form</b>	Liquid
<b>Clonality</b>	Polyclonal
<b>Alias</b>	BBC3, JFY1, JFY-1,
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	BBC3

## Image



Western blot analysis of extracts from 293 and MCF-7 cells using Puma Antibody.