





Phospho-SRC (Tyr418) Antibody

Product Code	CSB-PA014991
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12931
Immunogen	Peptide sequence around phosphorylation site of tyrosine 418(N-E-Y(p)-T-A)derived from Human Src.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of Src only when phosphorylated at tyrosine 418.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:100
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	ASV; SRC1; c-SRC; p60-Src
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SRC

Image



Western blot analysis of extracts from COLO205 cells using Src (Ab-418) antibody (Lane 1 and 2) and Src (phospho-Tyr418) antibody (Lane 3 and 4).



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Src (phospho- Tyr418) antibody.



Product Modify

Phospho-Tyr418

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.