




AFP Antibody

Product Code	CSB-PA075526
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02771
Immunogen	A synthesized peptide derived from Human Cyclin E1
Raised In	Rabbit
Species Reactivity	Human
Specificity	Cyclin E1 Antibody detects endogenous levels of total Cyclin E1
Tested Applications	ELISA,WB;WB:1:500-1:2000
Form	Supplied at 2mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.3, 0.05% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	alpha-fetoprotein, FETA, HPAFP
Immunogen Species	Homo sapiens (Human)
Target Names	AFP

Image



Western blot analysis of Cyclin E1 expression in Paclitaxel treated HeLa whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.