




CNTFR Antibody

Product Code	CSB-PA204475
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26992
Immunogen	Fusion protein corresponding to a region derived from 23-342 amino acids of Human ciliary neurotrophic factor receptor
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total CNTFR protein.
Tested Applications	ELISA,WB;WB:1:500-1:2000
Form	Supplied at 0.5mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.3, 0.05% sodium azide and 50% glycerol.
Clonality	Polyclonal
Alias	ciliary neurotrophic factor receptor
Immunogen Species	Homo sapiens (Human)
Target Names	CNTFR

Image



Gel: 10%SDS-PAGE Lysate: 40 ug Mouse skeletal muscle tissue lysate Primary antibody: 1/300 dilution Secondary antibody: Goat anti Rabbit IgG - H&L (HRP) at 1/10000 dilution Exposure time: 60 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.