




CLDN1 Antibody

Product Code	CSB-PA781428
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95832
Immunogen	A synthesized peptide derived from Human Claudin 1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	Claudin 1 antibody detects endogenous levels of total Claudin 1
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 °C.Stable for 12 months from date of receipt
Purification Method	Affinity-chromatography
Clonality	Polyclonal
Alias	claudin 1; Claudin-1; CLD1; CLDN1; ILVASC; SEMP1; Senescence-associated epithelial membrane protein; senescence-associated epithelial membrane protein 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLDN1

Image



Western blot analysis on HeLa cell lysate using Claudin 1 Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.