




MOK Antibody

Product Code	CSB-PA775498
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ07
Immunogen	A synthesized peptide derived from Human RAGE
Raised In	Rabbit
Species Reactivity	Human
Specificity	RAGE Antibody detects endogenous levels of RAGE
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:2000, IHC: 1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	MAPK/MAK/MRK overlapping kinase; MOK; MOK protein kinase; RAGE; RAGE-1; RAGE1; renal cell carcinoma antigen (MOK protein kinase); renal tumor antigen; Renal tumor antigen 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MOK

Image



Western blot analysis on HepG2 cell lysate using RAGE Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.