




TIMP1 Antibody

Product Code	CSB-PA020114
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01033
Immunogen	A synthesized peptide derived from Human TIMP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	TIMP1 Antibody detects endogenous levels of total TIMP1
Tested Applications	ELISA,WB;WB:1:500-1:2000
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 °C.Stable for 12 months from date of receipt
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	CLGI; Collagenase inhibitor; EPA; EPO; erythroid potentiating activity; Erythroid-potentiating activity; Fibroblast collagenase inhibitor; FLJ90373; HCl; Metalloproteinase inhibitor 1; TIMP; TIMP metallopeptidase inhibitor 1; TIMP-1; TIMP1; tissue inhibi
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TIMP1

Image



Western blot analysis of TIMP1 expression in 293 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.