




FOS Antibody

Product Code	CSB-PA138717
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01100
Immunogen	A synthesized peptide derived from Human Fos
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	Fos antibody detects endogenous levels of total Fos.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH7.4, 150mM NaCl, 0.02% sodium azide and50% glycerol.
Purification Method	Affinity-chromatography
Clonality	Polyclonal
Alias	activator protein 1; AP-1; C-FOS; cellular oncogene c-fos; Cellular oncogene fos; FBJ murine osteosarcoma viral (v-fos) oncogene homolog (oncogene FOS); FBJ murine osteosarcoma viral oncogene homolog; FOS; G0/G1 switch regulatory protein 7; G0S7; Proto-o
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOS

Image



Western blot analysis on K562 cell lysate using Fos antibody, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.