




# PLCG1 Antibody

<b>Product Code</b>	CSB-PA889574
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P19174
<b>Immunogen</b>	A synthesized peptide derived from Human PLCG1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	PLCG1 Antibody detects endogenous levels of total PLCG1
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
<b>Form</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 °C.Stable for 12 months from date of receipt
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Polyclonal
<b>Alias</b>	1-phosphatidylinositol-4,5-bisphosphate phosphodiesterase gamma 1, EC 3.1.4.11, PLC-148, PLC-II, PLC-gamma1, PLC1, PLCG-1, Phospholipase C-gamma-1
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PLCG1

## Image



Western blot analysis of PLCG1 expression in EGF treated 293 whole cell lysates. The lane on the left is treated with the antigen-specific peptide.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.