




PLCG2 Antibody

Product Code	CSB-PA890071
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16885
Immunogen	A synthesized peptide derived from Human PLCG2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	PLCG2 Antibody detects endogenous levels of total PLCG2
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	1-phosphatidylinositol-4,5-bisphosphate phosphodiesterase gamma 2, EC 3.1.4.11, PLC-IV, PLC-gamma2, Phospholipase C-gamma-2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLCG2

Image



Western blot analysis of PLCG2 expression in COLO205 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.