




PRKCD Antibody

Product Code	CSB-PA193274
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05655
Immunogen	A synthesized peptide derived from Human PKC delta
Raised In	Rabbit
Species Reactivity	Human
Specificity	PKC delta Antibody detects endogenous levels of total PKC delta
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	EC 2.7.11.13, KPCD, PKC-delta, PRKCD, kinase PKC-delta, nPKC-delta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCD

Image



Western blot analysis of PKC delta expression in MCF7 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.