




PRKCZ Antibody

Product Code	CSB-PA058972
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05513
Immunogen	A synthesized peptide derived from Human PKC zeta
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	PKC zeta Antibody detects endogenous levels of total PKC zeta
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	EC 2.7.11.13, KPCZ, PKC-zeta, PKC2, PRKCZ, Protein kinase C zeta type, kinase PKC-zeta, nPKC-zeta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCZ

Image



Western blot analysis of PKC zeta expression in PMA treated COS7 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.