




PRKCZ Antibody

Product Code	CSB-PA140454
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05513
Immunogen	A synthesized peptide derived from Human PKC theta
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	PKC theta Antibody detects endogenous levels of total PKC theta
Tested Applications	ELISA,WB;WB:1:500-1:2000
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	EC 2.7.11.13, KPCT, PKC-theta, PKCQ, PRKCQ, PRKCT, Protein kinase C, theta type, kinase PKC-theta, nPKC-theta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCZ

Image



Western blot analysis of PKC theta expression in PMA treated Jurkat whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.