




ELK1 Antibody

Product Code	CSB-PA255046
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19419
Immunogen	A synthesized peptide derived from Human Elk1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	Elk1 Antibody detects endogenous levels of total Elk1
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	ELK1, ETS-domain protein ELK-1, ETS-domain protein Elk-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELK1

Image



Western blot analysis of Elk1 expression in heatshock treated HeLa whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.