




TGFBR2 Antibody

Product Code	CSB-PA120767
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37173
Immunogen	A synthesized peptide derived from Human TGFR2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	TGFR2 antibody detects endogenous levels of total TGFR2
Tested Applications	ELISA,WB;WB:1:500-1:2000
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Clonality	Polyclonal
Alias	AAT3; FAA3; LDS1B; LDS2B; MFS2; RIIC; TAAD2; TbetaR-II; TGF-beta receptor type II; TGF-beta receptor type IIB; TGF-beta receptor type-2; TGF-beta type II receptor; TGFbeta-RII; TGFBR2; TGFR-2; TGFR2; transforming growth factor beta receptor type IIC; tra
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGFBR2

Image



Western blot analysis on NIH-3T3 cell lysate using TGFR2 Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.