




CREBBP Antibody

Product Code	CSB-PA027631
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92793
Immunogen	A synthesized peptide derived from Human CBP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	CBP Antibody detects endogenous levels of CBP
Tested Applications	ELISA, WB, IHC, IF; WB:1:500-1:2000, IHC:1:50-1:200, IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	CBP; CREB binding protein; CREB-binding protein; CREBBP; KAT3A; RSTS
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CREBBP

Image



Western blot analysis on HeLa cell lysate using CBP Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.