




Cleaved-CASP3 (A175) Antibody

Product Code	CSB-PA148736
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42574
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	Caspase 3 (p17,Cleaved-Asp175) Antibody detects endogenous levels of fragment of activated Caspase 3 resulting from cleavage adjacent to Asp175.
Tested Applications	ELISA, WB, IHC; WB:1:500-1:2000, IHC:1:50-1:200
Form	Liquid
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	Apopain precursor; CASP-3; CPP-32; CPP32; Caspase-3; Cysteine protease CPP32; ICE3; SCA-1; SREBP cleavage activity 1; Yama protein,p17
Immunogen Species	Homo sapiens (Human)
Target Names	CASP3

Image



Western blot analysis of extracts from COLO cells, using Caspase 3 (p17,Cleaved-Asp175) Antibody. The lane on the left is treated with the synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.