




# Cleaved-CASP3 (A175) Antibody

<b>Product Code</b>	CSB-PA148736
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P42574
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Specificity</b>	Caspase 3 (p17,Cleaved-Asp175) Antibody detects endogenous levels of fragment of activated Caspase 3 resulting from cleavage adjacent to Asp175.
<b>Tested Applications</b>	ELISA, WB, IHC; WB:1:500-1:2000, IHC:1:50-1:200
<b>Form</b>	Liquid
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Apopain precursor; CASP-3; CPP-32; CPP32; Caspase-3; Cysteine protease CPP32; ICE3; SCA-1; SREBP cleavage activity 1; Yama protein,p17
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CASP3

## Image



Western blot analysis of extracts from COLO cells, using Caspase 3 (p17,Cleaved-Asp175) Antibody. The lane on the left is treated with the synthesized peptide.