




ATP1A1 Antibody

Product Code	CSB-PA947783
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05023
Immunogen	A synthesized peptide derived from Human ATPase
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	ATPase Antibody detects endogenous levels of total ATPase
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	A1A1, AT1A1, ATP1A1, Na ⁺ /K ⁺ transporting, alpha 1 polypeptide, EC 3.6.3.9, Hypothetical protein MGC38419, Na ⁺ /K ⁺ ATPase 1, Na,K-ATPase 1, Sodium pump 1, Sodium/potassium-transporting ATPase alpha-1 chain precursor
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATP1A1

Image



Western blot analysis of ATPase expression in PMA treated 293 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.