




MAPK3/MAPK1 Antibody

Product Code	CSB-PA226662
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Swiss-Prot#: P27361/P28482
Immunogen	A synthesized peptide derived from Human ERK1/2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	ERK1/2 antibody detects endogenous levels of total ERK1/2
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Clonality	Polyclonal
Alias	Extracellular signal-regulated kinase 1, ERK-1, Insulin-stimulated MAP2 kinase, MAP kinase 1, MAPK 1, p44-ERK1, ERT2, p44-MAPK, Microtubule-associated protein 2 kinase extracellular signal-regulated kinase 1, ERK-1, Insulin-stimulated MAP2 kinase, MAP kin
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK3/MAPK1

Image



Western blot analysis on COLO205 cell lysate using ERK1/2 Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.