




HIF1A Antibody

Product Code	CSB-PA930291
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16665
Immunogen	A synthesized peptide derived from Human HIF1a.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	HIF1a antibody detects endogenous levels of HIF1a.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antigen Affinity Purified
Clonality	Polyclonal
Alias	HIF-1a,HIF1A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HIF1A

Image



Western blot analysis on LOVO using HIF1a antibody ,The lane on the left is blocked with the antigen-specific peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.