




HIF1A Antibody

Product Code	CSB-PA952643
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16665
Immunogen	Peptide sequence around aa.817-821(E-L-L-R-A) derived from Human HIF1A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total HIF1A protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	HIF1; MOP1; PASD8; HIF-1A; bHLHe78
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HIF1A

Image



Western blot analysis of extracts from K562 and HL-60 cells using HIF1A Antibody .

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.