




Phospho-INSR (Tyr1361) Antibody

Product Code	CSB-PA202534
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06213
Immunogen	Peptide sequence around phosphorylation site of Tyrosine 1361(I-P-Y(p)-T-H) derived from Human IR.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of IR only when phosphorylated at Tyrosine 1361.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	HHF5, CD220
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	INSR

Image



Western blot analysis of extracts from 293 cells treated with Heatstock using IR (Phospho-Tyr1361) Antibody. The lane on the right is treated with the antigen-specific peptide.

Product Modify	Phospho-Tyr1361
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.