




Phospho-MED1 (Thr1457) Antibody

Product Code	CSB-PA045259
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15648
Immunogen	Peptide sequence around phosphorylation site of Threonine 1457(A-Y-T(p)-P-Q) derived from Human PPAR-BP.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of PPAR-BP only when phosphorylated at Threonine 1457.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	PBP, CRSP1, RB18A, TRIP2, PPARBP
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MED1

Image



Western blot analysis of extracts from HUVEC cells treated with serum using PPAR-BP (Phospho-Thr1457) Antibody. The lane on the left is treated with the antigen-specific peptide.

Product Modify	Phospho-Thr1457
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.