




Phospho-ATRIP (Ser68) Antibody

Product Code	CSB-PA924373
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXE1
Immunogen	Peptide sequence around phosphorylation site of Serine 68(L-A-S(p)-Q-A) derived from Human ATRIP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of ATRIP only when phosphorylated at Serine 68.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Liquid
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	ATRIP
Immunogen Species	Homo sapiens (Human)
Target Names	ATRIP

Image



Western blot analysis of extracts from 3T3 cells using ATRIP (Phospho-Ser68) Antibody. The lane on the left is treated with the antigen-specific peptide.

Product Modify	Phospho-Ser68
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.