




Phospho-ERBB4 (Tyr1284) Antibody

Product Code	CSB-PA696402
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15303
Immunogen	Peptide sequence around phosphorylation site of Tyrosine 1284(P-E-Y(p)-L-S) derived from Human HER4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of p38 HER4 only when phosphorylated at Tyrosine 1284.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	HER4, p180erbB4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERBB4

Image



Western blot analysis of extracts from HUVEC cells treated with EGF using HER4 (Phospho-Tyr1284) Antibody. The lane on the left is treated with the antigen-specific peptide.

Product Modify Phospho-Tyr1284