




Phospho-FOXO4 (Thr451) Antibody

Product Code	CSB-PA939536
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98177
Immunogen	Peptide sequence around phosphorylation site of Threonine 451(L-G-T(p)-P-V) derived from Human FOXO4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of FOXO4 only when phosphorylated at Threonine 451.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using
Clonality	Polyclonal
Alias	AFX, AFX1, MLLT7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO4

Image



Western blot analysis of extracts from HUVEC cells treated with EGF using FOXO4 (Phospho-Thr451) Antibody. The lane on the left is treated with the antigen-specific peptide.

Product Modify	Phospho-Thr451
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.