




Phospho-ROCK2 (Tyr722) Antibody

Product Code	CSB-PA035451
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75116
Immunogen	Peptide sequence around phosphorylation site of Tyrosine 722(K-I-Y(p)-E-S) derived from Human ROCK2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of ROCK2 only when phosphorylated at Tyrosine 722.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	ROCK-II
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ROCK2

Image



Western blot analysis of extracts from Rat kidney using ROCK2 (Phospho-Tyr722) antibody. The lane on the left is treated with the antigen-specific peptide.

Product Modify	Phospho-Tyr722
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.