





ANGPTL3 Antibody

Product Code	CSB-PA169214
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5C1
Immunogen	Fusion protein corresponding to C terminal 300 amino acids of Human Angiopoietin-related protein 3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Liquid
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANGPTL3

Image



Gel: 10%SDS-PAGE Lysate: 30 ug Human liver cancer tissue lysate Primary antibody: 1/400 dilution Secondary antibody: Donkey anti Rabbit IgG - H&L (HRP) at 1/ 3000 dilution Exposure time: 1 minute



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using 31029(ANGPTL3 Antibody) at dilution 1/8, on the right is treated with the fusion protein.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.