




Phospho-GATA3 (Ser308) Antibody

Product Code	CSB-PA172946
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23771
Immunogen	Peptide sequence around phosphorylation site of Serine 308(R-L-S(p)-A-A) derived from Human GATA3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of GATA3 only when phosphorylated at Serine 308.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	HDR, HDRS
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GATA3

Image



Western blot analysis of extracts from HUVEC cells using GATA3 (Phospho-Ser308) Antibody. The lane on the right is treated with the antigen-specific peptide.

Product Modify	Phospho-Ser308
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.